



IONIZATION Series

PARTS LIST FOR OPERATION MAINTENANCE

<u>Part Number</u>	<u>Description</u>
410027	H25 Halogen Sensor
410018	Fan Filter
410159	28VDC Bulb (System 2000 & 3000 only)
410120	Pressure Gage
410136	Orifice .010
410144	Coarse Filter
410162	5 Micron Filter Replacement
410224	Diaphragm Pump (System 1000 only – 1030 and newer)
410107	GAST Pump (System 2000 & 3000 only)
420107	Service Kit for GAST Pump
410353	GAST Pump Filters (4 per pump)
410409	Carbon Vanes For GAST Pump (4 per pump)

COMMONLY ORDERED PARTS FOR TWO (2) YEARS OF MAINTENANCE

- 410162 5 Micron Filter Replacement – Filter is normally white. If the filter turns gray it needs changed. Recommend a change of this filter every six (6) months or if needed.
- 410144 Coarse Filter – Coarse filters (brass) are located at the end of each sampling point (zone). Coarse Filters should be checked and changed on a quarterly basis depending on environmental conditions.